## KFM Spinning Disk Confocal Microscope

3D Experience on your Desktop



Ξ



See the difference!

## **KFM Confocal Microskope Adaptor**

The ideal sensor for rapid acquisition of microscopic samples providing superior and detailed reproduction of all surface details within seconds.

- Microscopic technology with true depth discrimination provides superior surface acquisition right to the physical limit.
- Surface based video realtime acquisition principle for reliable results within seconds.
- Tough, maintenance free, long-lived construction.
- Extensive assortment of selected premium class objectives available.
- Ideal for measurement of topography,

KFM Confocal Mikroscope

geometry, form and position in 3D.

Principle of Operation

**Electronics** Configuration

Microscope

Light Source

Measurement Range

Camera

Objective

Video real time based confocal detection using a spinning nipkow disk. Tiny pin holes on the disk provide true depth discrimination, imaging only surface details within the focus of the objective. Surface parts outside the focus remain dark and do not interfere with the imaging and evaluation process. A dedicated authentic peak algorithm determines both the profile value and an omnifocal image.

Destinated especially for applications

in optical system engineering, materials science,

mechanical engineering and metrology

in semiconductor industry

health care

Integrated sensor and piezo controller, sockets for power supply, USB2.0 and RS232.

Internal bright field reflected light illumination, motorised switching into confocal mode. camera resolution 768x582 pixel<sup>2</sup>, up to 48 fps.

high power LED (505nm, cyan) illumination.

up to  $400\mu$ m (piezo axis), several mm (motorised Z axis) depending on microscope.

1 Ox	20x	40x	50x	100x
0.5	0.75	0,95	0.8	0.9
1.0	1.0	0.14	1.0	1.0
1210x908	605x445	302x223	243x181	121x91
1.58	0.8	0.40	0.32/0.42	0.158/0.37
0.01	0.003	0.002	0.002	0.001
	0.0			

aprox. 5kg ca 300mm x 140mm x 100mm

USB 2.0, RS232, Driver-DLL Complete integration into measurement and evaluation software Inspector piezo axis or motorised Z axis motorised X/Y axes image stitching software



Photo diode cell in a CMOS array 50 x 40 μm



Mico-lens arrav 125 x 125 µm



Chromium spheres on a roller 300 x 230 µm



tooth surface 40 x 30 µm



OPM Nobelstraße 7 76275 Ettlingen

49(0)7243/529942 Fone 49(0)7243/524480 Fax info@opm-messtechnik.de www.opm-messtechnik.de

Numerical Aperture Stand Off (mm) 1 Image Field ( $\mu$ m x  $\mu$ m) Lateral resolution ( $\mu$ m) Axial resolution  $(\mu m)$ 

Weight Dimensions

System Interfaces Software Integration

geometrical/optical resolution. <sup>2</sup> megapixel camera available Specifications subject to change

Options